

PLAT

Research Grade Tenecteplase

Catalog No.	DHB88401A DHB88401B	Quantity:	100 µg 1.0 mg
Alternate Names:	Tissue-type plasminogen activator		
Description:	Tenecteplase is a recombinant and site directed mutant version of human tissue plasminogen activator (tPA, EC 3.4.21.68) and is a member of the serine protease.		
Concentration:	1.0 mg/ml		
UniProt ID:	P00750		
Source:	Mammalian cell line		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Formulation:	0.01M PBS, pH 7.4		
Purification:	Ion exchange chromatography		
Purity:	> 95% as determined by SDS-PAGE		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

