

## PLAT

### Research Grade Alteplase

<b>Catalog No.</b>	DHB88402A DHB88402B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Tissue-type plasminogen activator, tPA, PLAT		
<b>Description:</b>	Alteplase is a recombinant human tissue plasminogen activator with a C-terminal 6X His tag.		
<b>UniProt ID:</b>	P00750		
<b>Source:</b>	Mammalian cell line		
<b>Concentration:</b>	1.0 mg/ml		
<b>Tag:</b>	C-terminal 6x His		
<b>Endotoxin:</b>	< 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Nickel chromatography		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for 1 week or as aliquots for up to 1 year at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

