

## NGF

### Research Grade Anti-Human NGF/Beta-NGF (MEDI1912)

<b>Catalog No.</b>	DHB91706A DHB91706B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Beta-nerve growth factor, Beta-NGF		
<b>Description:</b>	Affinity purified Research Grade Anti-Human NGF/Beta-NGF (MEDI1912) is a human monoclonal antibody which targets beta-nerve growth factor		
<b>UniProt ID (target):</b>	P01138		
<b>Specificity:</b>	Human Beta-nerve growth factor		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Human IgG		
<b>Clone:</b>	MEDI1912		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Formulation:</b>	0.01M PBS, pH 7.4.		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Endotoxin:</b>	≤ 0.001 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

