

TNF

Research Grade Placulumab

Catalog No.	DHB94407A DHB94407B	Quantity:	100 µg 1.0 mg
Alternate Names:	ART621, CEP-37247, PN0621, CAS: 945781-29-3		
Description:	Research Grade Placulumab is a human monoclonal antibody which targets Tumor necrosis factor		
UniProt ID (target):	P01375		
Specificity:	Human Tumor necrosis factor		
Source:	XtenCHO cell line		
Isotype:	Human [V-kappa]2-Fc		
INN:	Placulumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.