

IFNA1; IFNA13

Research Grade Sifalimumab

Catalog No.	DHB94501A DHB94501B	Quantity:	100 µg 1.0 mg
Alternate Names:	MDX-1103, MEDI-545, CAS: 1006877-41-3		
Description:	Research Grade Sifalimumab is a human monoclonal antibody which targets Interferon alpha-1/13		
UniProt ID (target):	P01562		
Specificity:	Human Interferon alpha-1/13		
Source:	XtenCHO cell line		
Isotype:	Humanized IgG1-kappa		
INN:	Sifalimumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

