

IL2RA

Research Grade Inolimomab

Catalog No.	DHB95804A DHB95804B	Quantity:	100 µg 1.0 mg
Alternate Names:	CAS: 152981-31-2		
Description:	Research Grade Inolimomab is a mouse monoclonal antibody which Interleukin-2 receptor subunit alpha		
UniProt ID (target):	P01589		
Specificity:	Human Interleukin-2 receptor subunit alpha		
Source:	XtenCHO cell line		
Isotype:	Mouse IgG1-kappa		
INN:	Inolimomab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

