

CD4

Research Grade Tregalizumab

Catalog No.	DHB95903A DHB95903B	Quantity:	100 µg 1.0 mg
Alternate Names:	BT-061, hB-F5, hBF5(H37L_L4M), CAS: 1207446-68-1		
Description:	Research Grade Tregalizumab is a humanized monoclonal antibody which targets T-cell surface glycoprotein CD4		
UniProt ID (target):	P01730		
Specificity:	Humanized T-cell surface glycoprotein CD4		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
INN:	Tregalizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

