

APOC3

Research Grade STT-5058 mAb

Catalog No.	DHB99001A DHB99001B	Quantity:	100 µg 1.0 mg
Alternate Names:	ARGX-116		
Description:	Research grade STT-5058 is a humanized monoclonal antibody that targets human Apolipoprotein C-III		
UniProt ID (target):	P02656		
Specificity:	Human Apolipoprotein C-III		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

