

FGB

Research Grade Biciromab

Catalog No.	DHB99201A DHB99201B	Quantity:	100 µg 1.0 mg
Alternate Names:	Indium (111In) biciromab, T2G1, CAS: 138783-13-8		
Description:	Research Grade Biciromab is a mouse monoclonal antibody which targets human Fibrinogen beta chain		
UniProt ID (target):	P02675		
Target:	Human Fibrinogen beta chain		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Biciromab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

