

TFRC

Research Grade Anti-Human CD71/TFRC (JST-TFR09)

Catalog No.	DHC02009A DHC02009B	Quantity:	100 µg 1.0 mg
Alternate Names:	JST-TFR09, JST TFR09, JSTTFR09		
Description:	Research Grade Anti-Human CD71/TFRC (JST-TFR09) is a human monoclonal antibody, which targets human Transferrin receptor protein 1 (CD71).		
UniProt ID (target):	P02786		
Specificity:	Human Transferrin receptor protein 1 (TfR1, CD71)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
Clone:	JST-TFR09		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

