

ERBB2

Research Grade Coprelotamab

Catalog No.	DHC09611A DHC09611B	Quantity:	100 µg 1.0 mg
Alternate Names:	GB-221, CAS: 2348469-84-9		
Description:	Research Grade Coprelotamab is a humanized monoclonal antibody directed against human Receptor tyrosine-protein kinase erbB-2 (CD340, HER2)		
UniProt ID (target):	P04626		
Target:	Human Receptor tyrosine-protein kinase erbB-2		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Coprelotamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.