

IL6

Research Grade Elsilimomab

Catalog No.	DHC15806A DHC15806B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-6, B-cell stimulatory factor 2, BSF-2, CTL differentiation factor, CDF, Hybridoma growth factor, Interferon beta-2, IFN-beta-2		
Description:	Elsilimomab (clone B-E8) is a mouse monoclonal antibody that targets and blocks human IL-6		
UniProt ID (target):	P05231		
Specificity:	Human IL-6		
Source:	XtenCHO cell line		
Isotype:	Mouse IgG1 kappa		
INN:	Elsilimomab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

