

## INSR

### Research Grade Valanafusp Alfa

<b>Catalog No.</b>	DHC19101A DHC19101B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	AGT-181, HIRMAb-IDUA, CAS: 1815583-32-4		
<b>Description:</b>	Research Grade Valanafusp Alfa is a chimeric monoclonal antibody directed against human insulin receptor (CD220)		
<b>UniProt ID:</b>	P06213		
<b>Target:</b>	Human Insulin receptor (CD220)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Chimeric IgG1 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>INN:</b>	Valanafusp Alfa		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or as aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

