

INSR

Research Grade Anti-Human CD220/INSR (IRAB-B)

Catalog No.	DHC19103A DHC19103B	Quantity:	100 µg 1.0 mg
Alternate Names:	IRAB-B, IRAB B, IRABB		
Description:	Research Grade Anti-Human CD220/INSR (IRAB-B) is an affinity purified human monoclonal antibody directed against human insulin receptor (CD220)		
UniProt ID (target):	P06213		
Target:	Human Insulin receptor (CD220)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G purified from cell culture supernatant		
Clone:	IRAB-B		
Application:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or as aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

