

CSF1R

Research Grade Cabiralizumab

Catalog No.	DHC25201A DHC25201B	Quantity:	100 µg 1.0 mg
Alternate Names:	FPA-008, CAS: 1613144-80-1		
Description:	Research Grade Cabiralizumab is a humanized monoclonal antibody which targets Macrophage colony-stimulating factor 1 receptor		
UniProt ID (target):	P07333		
Specificity:	Human Macrophage colony-stimulating factor 1 receptor		
Source:	XtenCHO cell line		
Isotype:	Human IgG4-kappa		
INN:	Cabiralizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

