

CD3E

Research Grade Otelixizumab

Catalog No.	DHC27702A DHC27702B	Quantity:	100 µg 1.0 mg
Alternate Names:	ChAglyCD, TRX4, CAS: 881191-44-2		
Description:	Research Grade Otelixizumab is a chimeric humanized monoclonal antibody which targets T-cell surface glycoprotein CD3 epsilon chain		
UniProt ID (target):	P07766		
Specificity:	Human T-cell surface glycoprotein CD3 epsilon chain		
Source:	XtenCHO cell line		
Isotype:	Chimeric Human IgG1-lambda		
INN:	Otelixizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

