

CD3E

Research Grade Dafsolimab

Catalog No.	DHC27707A DHC27707B	Quantity:	100 µg 1.0 mg
Alternate Names:	SPV-T3a, CAS: 2367001-71-4		
Description:	Research Grade Dafsolimab is a murine monoclonal antibody directed against human T-cell surface glycoprotein CD3 epsilon chain		
UniProt ID (target):	P07766		
Target:	Human T-cell surface glycoprotein CD3 epsilon chain		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgG2b kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Dafsolimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

