

IGF1R

Research Grade Dalotuzumab

Catalog No.	DHC29901A DHC29901B	Quantity:	100 µg 1.0 mg
Alternate Names:	MK-0646, CAS: 1005389-60-5		
Description:	Research Grade Dalotuzumab is a humanized monoclonal antibody which targets Human Insulin-like growth factor 1 receptor (CD221)		
UniProt ID (target):	P08069		
Specificity:	Human Insulin-like growth factor 1 receptor (CD221)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Dalotuzumab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

