

SERPINF2

Research Grade Anti-Human SERPINF2/ α 2-antiplasmin (DS9231) mAb

Catalog No.	DHC35601A DHC35601B	Quantity:	100 μ g 1.0 mg
Alternate Names:	Alpha-2-antiplasmin, Alpha-2-AP, Alpha-2-plasmin inhibitor, Alpha-2-PI, Serpin F2		
Description:	Research Grade Anti-Human SERPINF2/ α 2-antiplasmin (DS9231) is an affinity purified human monoclonal antibody directed against human alpha-2 antiplasmin.		
UniProt ID (target):	P08697		
Target:	Human alpha-2 antiplasmin		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 lambda		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	\leq 0.01 EU/ μ g by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Clone:	DS9231		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com