

CSF1

Research Grade Anti-Human CSF1/M-CSF (PD-0360324)

Catalog No.	DHC42103A DHC42103B	Quantity:	100 µg 1.0 mg
Alternate Names:	PD-0360324, PD 0360324, PD0360324		
Description:	Research Grade Anti-Human CSF1/M-CSF (PD-0360324) is an affinity purified human monoclonal antibody which targets Macrophage colony-stimulating factor 1		
UniProt ID (target):	P09603		
Specificity:	Human Macrophage colony-stimulating factor 1 (CSF-1)		
Source:	Mammalian		
Isotype:	Human IgG		
Clone:	PD-0360324		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Application:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

