


KIT

Research Grade Briquilimab

Catalog No.	DHC83102A DHC83102B	Quantity:	100 µg 1.0 mg
Alternate Names:	AMG 191, AMG-191, AMG191, JSP191, CAS: 2574591-89-0		
Description:	Research Grade Briquilimab is a chimeric humanized monoclonal antibody which targets human Mast/stem cell growth factor receptor Kit		
UniProt ID:	P10721		
Specificity:	Human Mast/stem cell growth factor receptor Kit (CD117)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Humanized IgG1 kappa		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.001 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Functional studies		
INN:	Briquilimab		
Manufacturer:	Antibody System		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

SDS-PAGE



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.