

FGFR1 & KLB Research Grade Fazpilodemab

Catalog No.	DHC88101A DHC88101B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, BFKB8488A, CAS: 2517935-02-1		
Description:	Research Grade Fazpilodemab is a bispecific chimeric humanized monoclonal antibody which targets human Fibroblast growth factor receptor 1 (CD331) and Beta-klotho		
UniProt ID:	P11362 & Q86Z14		
Specificity:	Human FGFR-1 (CD331) and Beta-klotho		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Fazpilodemab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.