

MS4A1

Research Grade Ocaratuzumab

Catalog No.	DHC90712A DHC90712B	Quantity:	100 µg 1.0 mg
Alternate Names:	AME-133v, LY2469298, CAS: 1169956-08-4		
Description:	Research Grade Ocaratuzumab is A humanized monoclonal antibody which targets B-lymphocyte antigen CD20		
UniProt ID (target):	P11836		
Specificity:	Human B-lymphocyte antigen CD20		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
INN:	Ocaratuzumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography ≤ 0.01 EU/µg by LAL analysis		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

