

MS4A1 & CD3E

Research Grade Mosunetuzumab

Catalog No.	DHC90716A DHC90716B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, BTCT4465A, RG-7828, RO7030816, CAS: 1905409-39-3		
Description:	Research Grade Mosunetuzumab is a human monoclonal antibody which targets B-lymphocyte antigen CD20 & T-cell surface glycoprotein CD3 epsilon chain		
UniProt ID (target):	P11836 & P07766		
Specificity:	Human B-lymphocyte antigen CD20 & T-cell surface glycoprotein CD3 epsilon chain		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
INN:	Mosunetuzumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

