

ITGA4 & ITGB7

Research Grade Etrolizumab

Catalog No.	DHC98803A DHC98803B	Quantity:	100 µg 1.0 mg
Alternate Names:	Anti-Beta7, PRO145223, RG7413, rhuMAb Beta7, CAS: 1044758-60-2		
Description:	Research Grade Etrolizumab is a humanized monoclonal antibody which targets Integrin alpha-4 & Integrin beta-7		
UniProt ID (target):	P13612 & P26010		
Specificity:	Human Integrin alpha-4 & Integrin beta-7		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
INN:	Etrolizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com