

CTLA4

Research Grade Ipilimumab

Catalog No.	DHD17201A DHD17201B	Quantity:	100 µg 1.0 mg
Alternate Names:	BMS-734016, MDX-010, MDX-CTLA-4, CAS: 477202-00-9		
Description:	Research Grade Ipilimumab is a human monoclonal antibody directed against Human Cytotoxic T-lymphocyte protein 4 (CTLA-4, CD152)		
UniProt ID (target):	P16410		
Target:	Human CTLA-4 (CD152)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Ipilimumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

