

## CTLA4

### Research Grade Tremelimumab

<b>Catalog No.</b>	DHD17203A DHD17203B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Ticilimumab, CP-675, CP-675, 206, CP-675206 clone 11.2.1, CAS: 745013-59-6		
<b>Description:</b>	Research Grade Tremelimumab is a human monoclonal antibody which targets Cytotoxic T-lymphocyte protein 4		
<b>UniProt ID (target):</b>	P16410		
<b>Specificity:</b>	Human Cytotoxic T-lymphocyte protein 4		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG2-kappa		
<b>INN:</b>	Tremelimumab		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Formulation:</b>	0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Endotoxin:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

