

CTLA4

Research Grade Nurulimab

Catalog No.	DHD17205A DHD17205B	Quantity:	100 µg 1.0 mg
Alternate Names:	BCD-145, CAS: 2168561-20-2		
Description:	Research Grade Nurulimab is a human monoclonal antibody which targets Cytotoxic T-lymphocyte protein 4		
UniProt ID (target):	P16410		
Specificity:	Human Cytotoxic T-lymphocyte protein 4		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
INN:	Nurulimab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

