

EPCAM

Research Grade Edrecolomab

| | | | |
|---------------------------------|---|------------------|------------------|
| Catalog No. | DHD17403A DHD17403B | Quantity: | 100 µg 1.0 mg |
| Alternate Names: | 17-1A, CAS: 156586-89-9 | | |
| Description: | Research Grade Edrecolomab is a mouse monoclonal antibody which targets human Epithelial cell adhesion molecule (CD326) | | |
| UniProt ID (target): | P16422 | | |
| Target: | Human Epithelial cell adhesion molecule (CD326) | | |
| Concentration: | 1.0 mg/ml, lot specific | | |
| Source: | XtenCHO | | |
| Isotype: | Mouse IgG2a kappa | | |
| Formulation: | 0.01M PBS buffer, pH 7.4 | | |
| Purity: | > 95% by SDS-PAGE | | |
| Endotoxin Level: | ≤ 0.01 EU/µg by LAL analysis | | |
| Purification: | Protein A/G affinity chromatography | | |
| INN: | Edrecolomab | | |
| Applications: | Functional studies | | |
| Storage & Stability: | Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

