

EPCAM

Research Grade Oportuzumab

Catalog No.	DHD17405A DHD17405B	Quantity:	100 µg 1.0 mg
Description:	Research Grade Oportuzumab is a humanized monoclonal antibody which targets human Epithelial cell adhesion molecule (CD326)		
UniProt ID (target):	P16422		
Target:	Human Epithelial cell adhesion molecule (CD326)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human scFv-kappa-heavy		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Oportuzumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

