

ITGB6

Research Grade Anti-Human ITGB6 (SGN-B6A)

Catalog No.	DHD29901A DHD29901B	Quantity:	100 µg 1.0 mg
Alternate Names:	SGN-B6A, SGN B6A, SGNB6A		
Description:	Research Grade Anti-Human ITGB6 (SGN-B6A) is a human monoclonal antibody which targets human Integrin beta-6		
UniProt ID (target):	P18564		
Target:	Human Integrin beta-6		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Clone:	SGN-B6A		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.