

LAG3

Research Grade Relatlimab

Catalog No.	DHD30403A DHD30403B	Quantity:	100 µg 1.0 mg
Alternate Names:	BMS-986016, ONO-4482, CAS: 1673516-98-7		
Description:	Research Grade Relatlimab is a human monoclonal antibody which targets human Lymphocyte activation gene 3 protein		
UniProt ID (target):	P18627		
Specificity:	Human Lymphocyte activation gene 3 protein (CD223)		
Source:	XtenCHO cell line		
Isotype:	Human IgG4 kappa		
INN:	Relatlimab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

