

LAG3

Research Grade Anti-Human CD223/LAG3 (GSK2831781)

Catalog No.	DHD30416A DHD30416B	Quantity:	100 µg 1.0 mg
Alternate Names:	Lymphocyte activation gene 3 protein, LAG-3, CD223, GSK2831781		
Description:	Research Grade Anti-Human CD223/LAG3 (GSK2831781) is a humanized monoclonal antibody that targets LAG-3+ activated T cells expressed in inflamed tissues.		
UniProt ID (target):	P18627		
Target:	Human Lymphocyte activation gene 3 protein (CD223)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells, expressed with glycosylated Fc		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Clone:	GSK2831781		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

