

TNFRSF1A

Research Grade Anti-Human CD120a/TNFRSF1A/TNFR1 (GSK2862277)

Catalog No.	DHD34004A DHD34004B	Quantity:	100 µg 1.0 mg
Alternate Names:	Tumor necrosis factor receptor superfamily member 1A, Tumor necrosis factor receptor 1, TNF-R1, TNFR1, CD120a, TBPI, p55, p60, CD120a, TNFRSF1A, GSK 2862277, GSK2862277		
Description:	Research Grade Anti-Human CD120a/TNFRSF1A/TNFR1 (GSK1995057) is a monoclonal antibody that targets human TNF-R1.		
UniProt ID (target):	P19438		
Specificity:	Human Tumor necrosis factor receptor 1 (TNFR1)		
Source:	XtenCHO cell line		
Isotype:	VHH-hFc		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Applications:	Functional studies		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com