

CD22

Research Grade Epratuzumab

Catalog No.	DHD38503A DHD38503B	Quantity:	100 µg 1.0 mg
Alternate Names:	90Y-DOTA-hLL2, 90Y-epratuzumab tetraxetan, AMG-412, 90Y-hLL2, hLL2, CAS: 205923-57-5		
Description:	Research Grade Epratuzumab is a human monoclonal antibody which targets B-cell receptor CD22		
UniProt ID (target):	P20273		
Specificity:	Human B-cell receptor CD22		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-nd		
INN:	Epratuzumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

