

CD22

Research Grade Suciraslimab

Catalog No.	DHD38509A DHD38509B	Quantity:	100 µg 1.0 mg
Alternate Names:	CAS: 2415439-63-1		
Description:	Research Grade Suciraslimab is a human monoclonal antibody, which targets human B-cell receptor CD22 (Siglec-2)		
UniProt ID:	P20273		
Specificity:	Human B-cell receptor CD22 (Siglec-2)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Suciraslimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

