

IL11

Research Grade Anti-Human IL11 (X203) mAb

Catalog No.	DHD42201A DHD42201B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-11, IL-11, EnX203, Adipogenesis inhibitory factor, AGIF, Oprelvekin		
Description:	Human monoclonal antibody (clone X203) which targets human Interleukin-11		
UniProt ID (target):	P20809		
Source:	XtenCHO		
Specificity:	Human IL-11		
Isotype:	Human IgG1 kappa		
Clone:	X203		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg		
Concentration:	1.0 mg/ml, lot specific		
Binding Assay:	EC50 = 202.9 ng/ml for anti-human IL-11 (X203) binding with human IL-11		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.