

## IL11

### Research Grade Anti-Human IL11 (BI 765423) mAb

<b>Catalog No.</b>	DHD42202A DHD42202B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Interleukin-11, IL-11, Adipogenesis inhibitory factor, AGIF, Oprelvekin, BI-765423, BI 765423, BI765423		
<b>Description:</b>	Human monoclonal antibody (BI 765423) which targets human Interleukin-11		
<b>UniProt ID (target):</b>	P20809		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Human IL-11		
<b>Isotype:</b>	Human IgG		
<b>Clone:</b>	BI 765423		
<b>Purification:</b>	Protein A or G affinity chromatography		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Endotoxin:</b>	< 0.01 EU/µg		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

