

ERBB3

Research Grade Barecetamab

Catalog No.	DHD48406A DHD48406B	Quantity:	100 µg 1.0 mg
Alternate Names:	ISU-104, ISU104, CAS: 2275727-74-5		
Description:	Research Grade Barecetamab is a human monoclonal antibody directed against human Receptor tyrosine-protein kinase erbB-3 (HER3)		
UniProt ID (target):	P21860		
Target:	Human Receptor tyrosine-protein kinase erbB-3 (HER3)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 lambda		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Barecetamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.