

TGM2

Research Grade Zampilimab

Catalog No.	DHD49201A DHD49201B	Quantity:	100 µg 1.0 mg
Alternate Names:	UCB-7858, UCB7858, CAS: 2098280-42-1		
Description:	Research Grade Zampilimab is a chimeric humanized monoclonal antibody directed against human Protein-glutamine gamma-glutamyltransferase 2		
UniProt ID (target):	P21980		
Target:	Human Protein-glutamine gamma-glutamyltransferase 2		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Zampilimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.