

CDH3

Research Grade Anti-Human CDH3/P-cadherin (PCA062) mAb

Catalog No.	DHD50503A DHD50503B	Quantity:	100 µg 1.0 mg
Alternate Names:	Cadherin-3, Placental cadherin (P-cadherin), CDH3, PCA062		
Description:	Research Grade Anti-Human CDH3/P-cadherin (PCA062) is a human monoclonal antibody (mAb) that targets Cadherin-3, a cell-to-cell adhesion molecule overexpressed and associated with cancer invasion and metastatic dissemination in many cancer types.		
UniProt ID (target):	P22223		
Specificity:	Human Cadherin-3		
Source:	Mammalian cell line		
Isotype:	Human IgG1		
Clone:	PCA062		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	Sterile filtered 0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Applications:	Functional studies		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com