

TNFSF4

Research Grade Amlitelimab

Catalog No.	DHD58802A DHD58802B	Quantity:	100 µg 1.0 mg
Alternate Names:	KY-1005, KY1005, KY1005-2D10, KY-1005-2D10, CAS: 2378692-15-8		
Description:	Research Grade Amlitelimab is a humanized monoclonal antibody directed against human Tumor necrosis factor ligand superfamily member 4 (OX40L, CD252)		
UniProt ID (target):	P23510		
Target:	Human Tumor necrosis factor ligand superfamily member 4 (OX40L, CD252)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Amlitelimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.