

KLRC1

Research Grade Monalizumab

Catalog No.	DHD71901A DHD71901B	Quantity:	100 µg 1.0 mg
Alternate Names:	IPH-2201, NN-8765, NNC-0141-0000-0100, NNC-141-01000, monalizumab, anti-NKG2A, humZ270, CAS: 1228763-95-8		
Description:	Research Grade Monalizumab is a humanized monoclonal antibody which targets NKG2-A/NKG2-B type II integral membrane protein (CD159a)		
UniProt ID (target):	P26715		
Specificity:	Human NKG2-A/NKG2-B type II integral membrane protein (CD159a)		
Source:	XtenCHO cell line		
Isotype:	Human IgG4 kappa		
INN:	Monalizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

