

KLRK1

Research Grade Tesnatilimab

Catalog No.	DHD72101A DHD72101B	Quantity:	100 µg 1.0 mg
Alternate Names:	JNJ-64304500, NN8555, NNC142-0002, JNJ 4500, JNJ4500, CAS: 2242758-08-1		
Description:	Research Grade Tesnatilimab is a human monoclonal antibody which targets human NKG2-D type II integral membrane protein (CD314)		
UniProt ID (target):	P26718		
Target:	Human NKG2-D type II integral membrane protein (CD314)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Tesnatilimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

