

## IL3RA

### Research Grade Talacotuzumab

<b>Catalog No.</b>	DHD72501A DHD72501B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	CSL-362-AML, CSL362, JNJ-56022473, CAS: 1826831-79-1		
<b>Description:</b>	Research Grade Talacotuzumab is a humanized monoclonal antibody which targets Interleukin-3 receptor subunit alpha.		
<b>UniProt ID (target):</b>	P26951		
<b>Specificity:</b>	Human Interleukin-3 receptor subunit alpha		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG1 kappa		
<b>INN:</b>	Vibecotamab		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Formulation:</b>	0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Endotoxin:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

