

IL3RA & CD3E

Research Grade Flotetuzumab

Catalog No.	DHD72504A DHD72504B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, MGD-006, MGD006, RES234, S80880, CAS: 1664355-28-5		
Description:	Research Grade Flotetuzumab is a bispecific humanized monoclonal antibody which targets Interleukin-3 receptor subunit alpha & T-cell surface glycoprotein CD3 epsilon chain		
UniProt ID (target):	P26951 & P07766		
Specificity:	Human Interleukin-3 receptor subunit alpha & T-cell surface glycoprotein CD3 epsilon chain		
Source:	XtenCHO cell line		
Isotype:	Human VH-VL-VH'-VL'		
INN:	Flotetuzumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.