

CD38

Research Grade Isatuximab

Catalog No.	DHD80802A DHD80802B	Quantity:	100 µg 1.0 mg
Alternate Names:	SAR-650984, hu38SB19, CAS: 1461640-62-9		
Description:	Research Grade Isatuximab is a chimeric monoclonal antibody directed against Human ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 (CD38)		
UniProt ID (target):	P28907		
Target:	Human ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 (CD38)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Isatuximab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

