

## CD38 & CD3E

### Research Grade Anti-Human CD38/ADPRC 1 & CD3e Bispecific Antibody (GBR-1342)

<b>Catalog No.</b>	DHD80822A DHD80822B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Bispecific, GBR-1342, ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1, 2'-phospho-ADP-ribosyl cyclase, Cyclic ADP-ribose hydrolase 1, cADPR hydrolase 1, T10, CD38, T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e		
<b>Description:</b>	Research Grade Anti-Human CD38/ADPRC 1 & CD3E Bispecific Antibody (GBR-1342) is directed against Cyclic ADP-ribose hydrolase 1 (CD38) and T-cell surface glycoprotein CD3 epsilon chain (CD3e)		
<b>UniProt ID (target):</b>	P28907 & P07766		
<b>Target:</b>	Human CD38 and CD3e		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Human half-IgG1-kappa/(scFv-kappa -heavy)-h-CH2-CH3		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	GBR-1342		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)