

AXL

Research Grade Mipasetamab

Catalog No.	DHD89503A DHD89503B	Quantity:	100 µg 1.0 mg
Alternate Names:	1H12-GalNAcN3, BGB-601, CAS: 2361055-48-1		
Description:	Research Grade Mipasetamab is a humanized monoclonal antibody directed against human Tyrosine-protein kinase receptor UFO (AXL oncogene)		
UniProt ID (target):	P30530		
Target:	Human Tyrosine-protein kinase receptor UFO (AXL oncogene)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Mipasetamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.