

CD70

Research Grade Anti-Human CD70 (SGN-CD70A)

Catalog No.	DHE01906A DHE01906B	Quantity:	100 µg 1.0 mg
Alternate Names:	SGN-CD70A, SGN CD70A, SGNCD70A		
Description:	Research Grade Anti-Human CD70 (SGN-CD70A) is a humanized monoclonal antibody affinity purified from cell culture supernatant.		
UniProt ID (target):	P32970		
Specificity:	Human CD70 antigen (TNFSF7, CD27 ligand)		
Source:	Mammalian cells		
Isotype:	Human IgG1-kappa		
Clone:	SGN-CD70A		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.